

Tru-Secure™

Peripheral Catheter Securement




- Provides for stabilization of peripheral catheters
- Protects insertion site with clear bio-occlusive dressing
- Replaces tape and Tegaderm™
- Non irritating foam adhesive safely secures and is easy to remove



Tru-Secure™ Peripheral Catheter Securement Device

The Tru-Secure™ Peripheral Catheter Securement Device makes securing a peripheral IV catheter easy. The single sterile pouched securement device takes the place of Tegaderm™ and tape reducing the amount of supplies needed to secure the catheter. Tru-Secure features a bio-occlusive covering to protect the insertion site while allowing for visualization of the access site and surrounding skin. The hypo-allergenic adhesive provides effective securement; yet does not require alcohol for removal. Tru-Secure is easy to apply and easy to remove making it less traumatic to the skin surface.



Ordering Information

| Catalog Number | Description | Quantity per Case |
|----------------|-----------------------|-------------------|
| TS-1000 | Tru-Secure Peripheral | 100 |



930 Wanamaker Ave
Ontario, CA 91761-8151

Toll-Free (800) 780-0012
Phone (909) 390-9080
Fax (909) 650-5894

www.maximusmedical.com